# 

# **Instrument Set for Arthroscopy System**



**ACL PCL Drill Guide System** 

Part No. M527001000



**Obturator for Drill Guide** 

Part No. M528001000

Cannulated Drill 4.5mm

**Part No.** M529000045

Flower Tip Reamer

Part No. M505000007

**Tunnel Notcher** 

Part No. M508001000

**3.5 Cannulated Screw Driver** 

**Part No.** M512001000

Tendon Stripper Dia 7 mm Open

Part No. |567001000

Tendon Stripper Dia 7 mm Closed

**Part No.** J567002000

\_\_\_\_

**Tendon Stripper Dia 5 mm Open** 

**Part No.** J567003000

Tendon Stripper Dia 5 mm Closed

**Part No.** J567004000

2.3mm Needle Drill Bit

Part No. J566000000

**Caliberated Probe** 

**Part No.** J501001000

**PCL Protector** 

Part No. |502001000

PCL Rasp

Part No. | 503000000

**Guide Wire: 9 inch** 

 Part No.
 Sizes

 J504001009
 1.2mm

 J505001009
 2mm

PITKARS PITKARS PITKARS 3



Tunnel Plug 8, 9, 10, 11, 12

Part No.

J506000000



**Graft Preparation Workstation** 

Part No.

J507001000



Rasp 20 Up

Part No.

1565001000



Rasp 20 Down

Part No.

J565002000

### **Arthroscopy System for Shoulder**



## **CUFFGRIP-T SUTURE ANCHOR**

Part No.	Dia
J165281000	2.8mm
J166351000	3.5mm
J167501000	5.0mm

# **Instruments for Rotogrip**



#### #ROTOGRIP

Part No.	Dia	Length
J115001000	4.75	24.5mm
J115002000	5.50	24.5mm

Following instruments are Also available

Probe 5.4mm

Punch 4.75mm & 5.50mm

Curette Both Sides Cut

Double Loaded Suture Needle

**Double Loaded Suture Handler** 

Shoulder Tissue elevator 15° & 30°

### Note:

- # Shoulder systems are sold together as an assembly of implant, suture & inserter.
- Inserter is single use only.